



## NRRL Medium No. 11

### Synthetischer nährstoffärmer [synthetic low-nutrient] Agar (SNA) w/carnation leaf

Per Liter of medium:

KH <sub>2</sub> PO <sub>4</sub> .....	1.0 g
KNO <sub>3</sub> .....	1.0 g
MgSO <sub>4</sub> •7H <sub>2</sub> O .....	0.5 g
KCl .....	0.5 g
Glucose .....	0.2 g
Sucrose .....	0.2 g
Agar .....	15.0g

Sterilize by autoclaving the above ingredients at 121C for 15 minutes, dispense in petri dishes. Place a piece of sterile carnation leaf (ca. 1 cm square) on the surface of the agar in each dish after the medium solidifies.