



NRRL Medium No. 14

Glycerol Asparagine Agar (ISP4)

Per Liter of medium:

Glycerol	10.0 g
Asparagine	1.0 g
K ₂ HPO ₄	1.0 g
Trace salts solution.....	1.0 ml
Agar	20.0 g

Adjust pH to 7.2 prior to autoclaving.

Trace Salts Solution

FeSO ₄ •7H ₂ O.....	0.1 g
MnCl ₂ •4H ₂ O.....	0.1 g
ZnSO ₄ •7H ₂ O.....	0.1 g
Water	100 ml